

Compounds & Plaster Selection Chart

Product Name	Compound Type	Working Time	Size	Jointing			Mechanical Tools	Sanding	Patching	Bonding	Wet Area	Skim Coat	Special Plaster	X-Ray	
				1st Coat	2nd Coat	Finishing Coat									
BaseCote®	Setting	45, 60 & 90 mins	20kg Bag	✓	✓		✓								
RediBase™	Air Drying	–	18kg Pail	✓	✓		✓								
GIB X-BLOCK® Note: to be used in conjunction with X-Block® Plasterboard.	Air Drying	–	25kg Pail	✓	✓		✓								✓
LiteFinish™	Air Drying	–	18kg Pail			✓	✓	Very Easy							
FinalCote®	Air Drying	–	20kg Pail			✓	✓	Easy							
TopCote 550®	Air Drying	–	20kg Pail			✓	✓	Easy							
All Purpose Premix	Air Drying	–	18kg Pail	✓	✓	✓	✓	Moderate	✓						
PatchCote 20™	Setting	20 mins	8kg Pail					Moderate	✓	✓					
Cornice Adhesive	Setting	45 & 60 mins	20kg Bag					Moderate	✓	✓					
Masonry Adhesive	Setting	65 mins	20kg Bag							✓					
Back-Blocking Adhesive	Setting	120 mins	20kg Bag							✓					
Premium Bond Stud Adhesive	Air Drying	–	5kg Pail, 600ml Sausage							✓					
Wet Area Sealant	Air Drying	–	450g Cartridge								✓				
Wet Area SealCote	Air Drying	–	4lt Tub								✓				
Wet Area Taping Cement	Air Drying	–	3.2 & 6.4kg Pail	✓	✓						✓				
Hard Finish	Setting	20 mins 70 mins (WA) 45 mins (SA)	20kg Bag									✓			
Ceraplast™	Setting	20 mins	20kg Bag											✓	
Casting Plaster	Setting	25 mins	20kg Bag											✓	
Dental Plaster	Setting	5 mins	20kg Bag											✓	
Pottery Plaster	Setting	10 mins	20kg Bag											✓	